



CRARY SHOES
CUSTOM SHOE ORDERING

DUNS #037083706

*Registered on SAM.gov as:
Cragon Inc. DBA Crary Shoe Co.*

Crary Shoes

Custom Molded Shoes Since 1978

T: 503-253-8984 | F: 503-253-2094 | admin@craryshoes.com | www.craryshoes.com

Shoe Ordering Process

1. Patient evaluation with the ordering physician

- A. Casting (STS casting socks only) on a PLATFORM covered in foam
- B. Measurements taken over socks, brace, bandages, etc.
(Please use metric – NOT inches)
- C. Fill out order form completely
- D. Photos (optional) – can be emailed to admin@craryshoes.com

2. Mail to :

Crary Shoes - ATTN: CMS
12035 NE Glenn Widing Dr.
Portland, OR 97220

3. We will review order and call with any questions.

Once your order has been received with all of the necessary information, we will ship your order in approximately 4-5 weeks. More complicated orders can take longer. Orders can be rushed for a fee. Please note, if an order is incomplete or needs a prior authorization/purchase order confirmation, this will delay the delivery time.

Call or email Crary Shoes with questions
that you don't find answers to in this
packet. **503.253.8984** or
admin@craryshoes.com

Crary Shoes

Custom Molded Shoes Since 1978

T: 503-253-8984 | F: 503-253-2094 | admin@craryshoes.com | www.craryshoes.com

STS Casting Socks

Fitted Sock Application Instructions

1. With the foot at 90° to the leg, place cutting strip on dorsum of foot.



2. Tape Channel tube on top of cutting strip.



3. Place plastic baggie of appropriate length on foot and leg.



4. Dip sock in cool water. Then gather moistened STS Sock and slip on foot.



5. Completely apply the sock to foot and leg. Stretch out any wrinkles, leave no overlapping, thoroughly massage sock to



STS Casting Socks

Fitted Sock Application Instructions

Continued

6.Place foot in desired neutral position, (Fig. 6a) or in a semi-weight bearing foot position (Fig. 6b).



7.Mark malleoli, 1st and 5th metatarsals heads or boney prominences with indelible pen.



8.Remove channel tube.



9.Cut sock with scissors or letter opener overcutting strip. F



10.Spread cured STS sock and remove from foot and leg. shape of foot, ankle and leg.



STS Casting Socks

Sizing Information	ANKLE CASTING SOCK		
	Item Number	Size	Shoe Size
	900-S	Small	W: 4-6 1/2
	900-M	Medium	W: 7-10, M: 7-9
	900-L	Large	W: 10 1/2 - 12, M: 9 1/2 - 11 1/2
	900-XL	X-Large	M: 12-15

Marking the Cast

& Work Order for Modifications

You can use an indelible pencil, but we find it easy to simply use a sharpie pen.

1. A Series of Straight Lines

Making this marking to show where you want material removed, which will result in additional support. Usually these areas are under the medial longitudinal arch or just behind the metatarsal heads for built in metatarsal support.



2. A Circle with an X in the Middle

Mark a simple X at the apex of a bony prominence. This is commonly used for hammer-toes, a prominent styloid or metatarsal head, etc. Draw a circle around the area. A 6mm built up will be applied at the X and will be tapered down to the edges of the circle.

3. A Starburst with a Dot in the Center

This mark is for a severe prominence. This area will be built up 10mm from the apex and tapered down. This is usually at a prominent navicular, cuboid, severe hammer toe, or metatarsal head.



4. Mark AFO/Braces

Please use the pen to outline the edges of the brace. This is extremely helpful when we are building the last so we don't make the shoe too long, too short depending on whether the brace is 3/4 or full length or too narrow at the heel if there are hinges.

IMPORTANT FOR PLANTAR SURFACE BUILDUPS:

When the buildups are needed on the plantar surface of the foot, please also mark these during the weight bearing tracing on the paper. This helps us to see the elongation of the foot. This gives us better accuracy for the location of the buildups.

Line Sheet Women's



Annie

COLORS: *BONE ONLY*

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain leather upper and breathable cloth lining



Addison (t-strap)

COLORS: **Black only**

- Adjustable velcro strap with d-ring on medial and lateral side
- Padded collar for comfort
- Full grain leather upper and leather lining
- Can be ordered with our without the leather flower detail
- ****NOT RECOMMENDED FOR SEVERE DEFORMITIES**



Lo-Sneaks_VELCRO

COLORS: **Black only**

- Dual Velcro straps w/ d-ring
 - Padded collars and tongue for comfort
 - Full grain leather upper and breathable cloth lining
 - Perforated leather over vamp for breathability and style
- sole can be black or white-
black is the default



Athletic

COLORS: **Black only**

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain Italian leather upper and breathable mesh lining
- Perforated leather over vamp for additional breathability and style



Lo-Sneaks

COLORS: **White/Grey only**

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain Italian leather upper and breathable mesh lining
- Perforated leather over vamp for additional breathability and style



Hi-Sneaks _lace or velcro

COLORS: **Black only**

- Just over the ankle for additional support
- Lace closure system
- Padded collars and tongue for comfort
- Full grain leather upper and breathable cloth lining
-

Crary Shoes

Custom Molded Shoes Since 1978

T: 503-253-8984 | F: 503-253-2094 | admin@craryshoes.com | www.craryshoes.com

Line Sheet Men's



Athletic

COLORS: Black only

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain Italian leather upper and breathable mesh lining
- Perforated leather over vamp for additional breathability and style



Explorer

COLORS: Brown OR White

- Dual Velcro straps w/ d-ring
- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain leather upper and breathable cloth lining
-



Lo-Sneaks VELCRO

COLORS: Black only

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain Italian leather upper and breathable mesh lining
- Perforated leather over vamp for additional breathability and style



Lo-Sneaks

COLORS: White/Grey only

- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain Italian leather upper and breathable mesh lining
- Perforated leather over vamp for additional breathability and style



Hi-Sneaks

COLORS: Black only

- Just over the ankle for additional support
- Lace closure system
- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain leather upper and breathable cloth lining



Hi-Sneaks Velcro

COLORS: Black only

- Just over the ankle for additional support
- Elastic over the instep and a single Velcro closure around the ankle
- Deep throat for easy entry
- Padded collars and tongue for comfort
- Full grain leather upper and breathable cloth lining

Crary Shoes

Custom Molded Shoes Since 1978

T: 503-253-8984 | F: 503-253-2094 | admin@craryshoes.com | www.craryshoes.com